

Product Information

Recombinant Anti-Human NGF Antibody Fab Fragment

Cat. No.: **MOM-18110-F(P)**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Human Antibody Fab Fragment binds selectively to Human NGF, expressed in E. coli

Antigen Description

Nerve growth factor is important for the development and maintenance of the sympathetic and sensory nervous systems. It stimulates division and differentiation of sympathetic and embryonic sensory neurons.

Specific Activity

Tested positive against native antigen.

Target

NGF

Source

Human

Species Reactivity

Human

Type

Fab Fragment based on Human IgG2 - kappa

Expression Host

E. coli

Purity

>95.0%, determined by analysis by RP-HPLC & analysis by SDS-PAGE.

Applications

Suitable for use in FC, IP, ELISA, Neut, FuncS, IF and most other immunological methods.

Storage

Store it under sterile conditions at -20°C upon receiving. Recommend to pack the protein into smaller quantities for optimal storage.

ANTIGEN GENE INFORMATION

Gene Name

[NGF nerve growth factor \(beta polypeptide\) \[Homo sapiens \]](#)

Official Symbol

NGF

Synonyms

NGF; nerve growth factor (beta polypeptide); NGFB; beta-nerve growth factor; nerve growth factor, beta subunit; HSN5; Beta-NGF; MGC161426; MGC161428;

Gene ID

[4803](#)

mRNA Refseq

[NM_002506](#)

Protein Refseq

[NP_002497](#)

MIM

[162030](#)

UniProt ID

P01138

Chromosome Location

1p13.1

Pathway

ARMS-mediated activation, organism-specific biosystem; Activation of TRKA receptors, organism-specific biosystem; Apoptosis, organism-specific biosystem; Apoptosis, conserved biosystem; Axonal growth stimulation, organism-specific biosystem; Cell death signalling via NRAGE, NRIF and NADE, organism-specific biosystem; Ceramide signalling, organism-specific biosystem;

Function

growth factor activity; nerve growth factor receptor binding; receptor signaling protein activity;