

Product Information

Recombinant Anti-Human CD6 Antibody

Cat. No.: **MOM-18247**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Humanized (from mouse) Antibody is against Human CD6, expressed in Chinese Hamster Ovary cells(CHO)

Antigen Description

This gene encodes a protein found on the outer membrane of T-lymphocytes as well as some other immune cells. The encoded protein contains three scavenger receptor cysteine-rich (SRCR) domains and a binding site for an activated leukocyte cell adhesion molecule. The gene product is important for continuation of T cell activation. This gene may be associated with susceptibility to multiple sclerosis (PMID: 19525953, 21849685). Multiple transcript variants encoding different isoforms have been found for this gene.

Specific Activity

Tested positive against native antigen.

Target

CD6

Immunogen

The details of the immunogen for this antibody are not available.

Source

Humanized (from mouse)

Species Reactivity

Human

Type

Humanized (from mouse) IgG1 - kappa

Expression Host

CHO

Predicted N terminal

H chain: EVQLVES; L Chain: DIQMTQS

Purity

>97%, by SDS-PAGE under reducing conditions and visualized by silver stain.

Applications

Suitable for use in IF, IP, Neut, FuncS, ELISA, FC, WB and most other immunological methods.

Storage

Store it under sterile conditions at -20°C upon receiving. Recommend to pack the protein into smaller quantities for optimal storage.

ANTIGEN GENE INFORMATION

Gene Name

[CD6 CD6 molecule \[Homo sapiens \]](#)

Official Symbol

CD6

Synonyms

CD6; CD6 molecule; CD6 antigen; T-cell differentiation antigen CD6; Tp120; T12; TP120; FLJ44171;

Gene ID

[923](#)

mRNA Refseq

[NM_001254750](#)

Protein Refseq

[NP_001241679](#)

MIM

[186720](#)

UniProt ID

P30203

Chromosome Location

11q12.2

Pathway

Cell adhesion molecules (CAMs), organism-specific biosystem; Cell adhesion molecules (CAMs), conserved biosystem;

Function

scavenger receptor activity;