

Product Information

Recombinant Anti-Human msln Antibody scFv Fragment

Cat. No.: **MOM-18592-S(P)**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Mouse Antibody scFv Fragment is directed against Human MSLN, expressed in E. coli

Antigen Description

Membrane-anchored forms may play a role in cellular adhesion. Megakaryocyte-potentiating factor (MPF) potentiates megakaryocyte colony formation in vitro.

Specific Activity

Tested positive against native antigen.

Target

MSLN

Immunogen

The details of the immunogen for this antibody are not available.

Source

Mouse

Species Reactivity

Human

Type

scFv

Expression Host

E. coli

Purity

>95.0% as determined by analysis by RP-HPLC.

Applications

Suitable for use in ELISA, WB, Neut and most other immunological methods.

Storage

Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze/thaw cycles.

ANTIGEN GENE INFORMATION

Gene Name

[MSLN mesothelin \[Homo sapiens \]](#)

Official Symbol

MSLN

Synonyms

MSLN; mesothelin; CAK1; MPF; CAK1 antigen; megakaryocyte potentiating factor; soluble MPF mesothelin related protein; pre-pro-megakaryocyte-potentiating factor; SMRP;

Gene ID

[10232](#)

mRNA Refseq

[NM_001177355](#)

Protein Refseq

[NP_001170826](#)

MIM

[601051](#)

UniProt ID

Q13421

Chromosome Location

16p13.3