

Product Information

Recombinant Anti-Human tnfsf12 Antibody

Cat. No.: **MOM-18627**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Mouse Antibody is specific to Human TNFSF12, expressed in Chinese Hamster Ovary cells(CHO)

Antigen Description

Binds to FN14 and possibly also to TNFRSF12/APO3. Weak inducer of apoptosis in some cell types. Mediates NF-kappa-B activation. Promotes angiogenesis and the proliferation of endothelial cells. Also involved in induction of inflammatory cytokines.

Specific Activity

Tested positive against native antigen.

Target

TNFSF12

Immunogen

Recombinant full length TWEAK (Human).

Source

Mouse

Species Reactivity

Human

Type

IgG

Expression Host

CHO

Purity

>95%, by SDS-PAGE with silver staining, under reducing conditions.

Applications

Suitable for use in Neut, FuncS, IF, ICC and most other immunological methods.

Storage

Store at -20°C. Avoid multiple freeze/thaw cycles.

ANTIGEN GENE INFORMATION

Gene Name

[TNFSF12 tumor necrosis factor \(ligand\) superfamily, member 12 \[Homo sapiens \]](#)

Official Symbol

TNFSF12

Synonyms

TNFSF12; tumor necrosis factor (ligand) superfamily, member 12; tumor necrosis factor ligand superfamily member 12; APO3L; DR3LG; TWEAK; APO3 ligand; APO3/DR3 ligand; TNF-related WEAKE inducer of apoptosis; MGC20669; MGC129581;

Gene ID

[8742](#)

mRNA Refseq

[NM_003809](#)

Protein Refseq

[NP_003800](#)

MIM

[602695](#)

UniProt ID

O43508

Chromosome Location

17p13.1

Pathway

Cytokine-cytokine receptor interaction, organism-specific biosystem; Cytokine-cytokine receptor interaction, conserved biosystem; Gene Expression, organism-specific biosystem; Regulation of mRNA Stability by Proteins that Bind AU-rich Elements, organism-specific biosystem; Stabilization of mRNA by HuR, organism-specific biosystem;

Function

cytokine activity; cytokine activity; protein binding; receptor binding; tumor necrosis factor receptor binding;