

Product Information

Recombinant Anti-Human IL17A Antibody scFv Fragment

Cat. No.: **MOM-H57-S(P)**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Humanized (from mouse) Antibody scFv Fragment is bind to Human IL17A, expressed in E. coli

Antigen Description

The protein encoded by this gene is a proinflammatory cytokine produced by activated T cells. This cytokine regulates the activities of NF-kappaB and mitogen-activated protein kinases. This cytokine can stimulate the expression of IL6 and cyclooxygenase-2

Specific Activity

IL17A (interleukin 17A, IL-17A) [Homo sapiens]

Target

IL17A

Source

Humanized (from mouse)

Species Reactivity

Human

Type

Humanized (from mouse) scFv

Expression Host

E. coli

Purity

>95%, by SDS-PAGE with silver staining, under reducing conditions.

Purification

Purified by Nickel ion affinity chromatography

Applications

Suitable for use in FC, IP, ELISA, Neut, FuncS, IF and most other immunological methods.

Cellular Localization

kappa

Storage

Store the antibody (in aliquots) at -20°C. Avoid repeated freezing and thawing of samples.

ANTIGEN GENE INFORMATION

Gene Name

[IL17A interleukin 17A \[Homo sapiens \]](#)

Official Symbol

IL17A

Synonyms

IL17A; interleukin 17A; CTLA8, IL17, interleukin 17 (cytotoxic T lymphocyte associated serine esterase 8); interleukin-17A; cytotoxic T lymphocyte associated protein 8; IL 17; IL 17A; CTLA-8; cytotoxic T-lymphocyte-associated antigen 8; cytotoxic T-lymphocyte-associated protein 8; cytotoxic T-lymphocyte-associated serine esterase 8; interleukin 17 (cytotoxic T-lymphocyte-associated serine esterase 8); IL17; CTLA8; IL-17; IL-17A;

Gene ID

[3605](#)

mRNA Refseq

[NM_002190](#)

Protein Refseq

[NP_002181](#)

MIM

[603149](#)

UniProt ID

Q16552

Chromosome Location

6p12

Pathway

Cytokine-cytokine receptor interaction, organism-specific biosystem; Cytokine-cytokine receptor interaction, conserved biosystem; IL23-mediated signaling events, organism-specific biosystem; IL27-mediated signaling events, organism-specific biosystem; Rheumatoid arthritis, organism-specific biosystem; Rheumatoid arthritis, conserved biosystem;

Function

cytokine activity;