

Product Information

Recombinant Anti-Human ITGA4+ITGB7 Antibody

Cat. No.: **MOM-H75**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Humanized (from mouse) Antibody is bind to Human ITGA4+ITGB7, expressed in HEK293

Specific Activity

ITGA4_ITGB7(integrin alpha4 (CD49d)_beta7, integrin alpha4_beta7, integrin α 4 β 7, lymphocyte Peyer s patch adhesion molecule 1, LPAM-1) [Homo sapiens]

Target

ITGA4+ITGB7

Source

Humanized (from mouse)

Species Reactivity

Human

Type

Humanized (from mouse) IgG1 - kappa

Expression Host

HEK293

Purity

>95%, by SDS-PAGE with silver staining, under reducing conditions.

Purification

Protein A affinity purified from an animal origin - free culture supernatant

Applications

Suitable for use in ELISA, WB, Neut and most other immunological methods.

Cellular Localization

kappa

Storage

Store at -20°C for long-term storage. Store at 2-8°C for up to one month. Avoid freeze/thaw cycles.

BACKGROUND

Introduction

Integrins are transmembrane receptors that are the bridges for cell-cell and cell-extracellular matrix (ECM) interactions. When triggered, integrins in turn trigger chemical pathways to the interior (signal transduction), such as the chemical

composition and mechanical status of the ECM, which results in a response (activation of transcription) such as regulation of the cell cycle, cell shape, and/or motility; or new receptors being added to the cell membrane

Keywords

ITGA4; ITGB7; Integrin $\alpha 4\beta 7$; CD49d/ $\beta 7$