

Product Information

MemDX™ Antibody Discovery - Human IGF-I R / CD221 (31-932) Membrane Protein, Partial, -His tag

Cat. No.: **MP0246F**

This product is for research use only and is not intended for diagnostic use.

This membrane protein is Human IGF-I R / CD221 (31-932). It has been tested in SDS-PAGE, ELISA. We provide this protein to facilitate your membrane protein antibody discovery and development.

Product Specifications

Host Species

Human

Target Protein

IGF-I R / CD221

Protein Length

ECD

Molecular Weight

This protein contains a furin convertase cleavage site, 737-RKRR-740, and will be partially processed into N (α chain) and C-terminal fragment (partial β chain) with calculated MW of 81.0 kDa and 23.8 kDa respectively. The protein migrates as 36 kDa (partial β chain), 96-116 kDa (α chain) and 120-140 kDa (α chain & partial β chain) due to glycosylation.

Sequence

AA Glu 31 - Asn 932 (Accession # P08069-1).

Product Description

Activity

Yes

Application

SDS-PAGE, ELISA

Expression Systems

HEK293

Tag

His tag at the C-terminus

Protein Format

Soluble

Form

LYOPH

Reconstitution

Please see Certificate of Analysis for specific instructions.

Endotoxin

<1.0 EU/μg by the LAL method

Purity

>90% as determined by SDS-PAGE.

Buffer

Lyophilized from 0.22 μm filtered solution in PBS, pH7.4. Normally trehalose is added as protectant before lyophilization.

Storage

Stored at lyophilized form at -20°C or lower. Avoid repeated freeze-thaw cycles.

The antigen can be stable for 12 months in lyophilized form after storage at -20°C to -80°C, 3 months under sterile conditions after reconstitution after storage at -80°C.

Target**Target Protein**

IGF-I R / CD221

Full Name

insulin like growth factor 1 receptor

Introduction

This receptor binds insulin-like growth factor with a high affinity. It has tyrosine kinase activity. The insulin-like growth factor I receptor plays a critical role in transformation events. Cleavage of the precursor generates alpha and beta subunits. It is highly overexpressed in most malignant tissues where it functions as an anti-apoptotic agent by enhancing cell survival. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.

Alternative Names

IGFR; CD221; IGFI; JTK13; insulin-like growth factor 1 receptor; IGF-I receptor; soluble IGF1R variant 1; soluble IGF1R variant 2

Gene ID

[3480](#)

UniProt ID

[P08069](#)