

Product Information

MemDX™ Membrane Protein Human ADGRF5 (Adhesion G protein-coupled receptor F5) for Antibody Discovery

Cat. No.: **MP0439X**

This product is for research use only and is not intended for diagnostic use.

This product is a 148.06 kDa Human ADGRF5 membrane protein expressed in *in vitro* wheat germ expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

ADGRF5

Protein Length

Full-length

Molecular Weight

148.06 kDa

TMD

7

Sequence

MKSPRRRTLCLMFIVYSSKAALNWNYESITHPLSLHEHEPAGEEALRQKRAVATKSPTAEEYTVNIEISFENASFLDPIKAYLNSLSFF

Product Description

Application

Antibody Production

Expression Systems

in vitro wheat germ expression system

Tag

NO

Protein Format

Liposome

Form

Liquid

Purification

None

Buffer

25 mM Tris-HCl of pH8.0 containing 2% glycerol

Storage

Store at +4°C for up to one week or several months at -80°C

Target**Target Protein**

ADGRF5

Full Name

Adhesion G protein-coupled receptor F5

Introduction

The ADGRF5 gene is conserved in chimpanzee, Rhesus monkey, dog, cow, mouse, rat, chicken, and frog.

Alternative Names

GPR116; KPG_001

Gene ID

[221395](#)

UniProt ID

[Q8IZF2](#)