

Product Information

MemDX™ Membrane Protein Human GPR135 (G protein-coupled receptor 135) Full Length

Cat. No.: **MPC0120K**

This product is for research use only and is not intended for diagnostic use.

This product is a 51.7 kDa Human GPR135 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

GPR135

Protein Length

Full length

Protein Class

GPCR

Molecular Weight

51.7 kDa

TMD

7

Sequence

MEEPQPPRPPASMLLGSQHSGAPSAAGPPGGTSSAATAAVLSFSTVATA
ALGNLSDASGGGTAAAPGGGGLGGGAAREAGAAVRRPLGPEAAPLLSHG
AAVAAQALVLLIFLLSSLGNCVMGVIVKHRQLRTVTNAFILSLSLSDL
LTALLCLPAAFLDLFTPPGGSAPAAAAGPWRGFCAASRFFSSCFGIVSTL
SVALISLDTRYCAIVRPPREKIGRRRALQLLAGAWLTALGFSLPWELLGAP
RELAAAQSFHGCLYRTSPDPAQLGAAFSVGLVVACYLLPFLLMCFCHYHI
CKTVRLSDVRVRPVNTYARVLRFFSEVRTATTVLIMIVFVICCWGPYCFL
VLLAAARQAQTMQAPSLLSVVAVWLTWANGAINPVIYAIRNPINISMLLGR
NREEGYRTRNVDAFLPSQGPGLQARSRSRLRNRYANRLGACNRMSSSNPA
SGVAGDVAMWARKNPVVLFCREGPPEPVTAVTKQPKSEAGDTSL

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

GPR135

Full Name

G protein-coupled receptor 135

Alternative Names

HUMNP1Y20; G-protein coupled receptor 135; probable G-protein coupled receptor 135; GPR135; G protein-coupled receptor 135

Gene ID

[64582](#)

UniProt ID

[Q8IZ08](#)