

Product Information

MemDX™ Membrane Protein Human TAAR3P (Trace amine associated receptor 3, pseudogene) Full Length

Cat. No.: **MPC0249K**

This product is for research use only and is not intended for diagnostic use.

This product is a 39 kDa Human TAAR3P membrane protein expressed in Baculovirus/Insect expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TAAR3P

Protein Length

Full length

Protein Class

GPCR

Molecular Weight

39 kDa

TMD

7

Sequence

MDLTYIPEDLSSCPKFVNKILSSHQPLFSCPGDNVFGYDWSHDYPLFGNL
VIMVSISHFKQLHSPTNFLILSMATTDFLLG FVIMPYSIMRSVESCWYFG
DGFCKFHTSFDMMLRLTSIFHLCSIAIDRFYAVCYPLHYTTKMTNSTIKQ
LLAFCWSVPALFSFGLVLSEADVSGMQSYKILVACFNFCALTFNKFVGTI
LFTTCFFTPGSIMVGIYGKIFIVSKQHARVISHVPENTKGAVKKHLSKKK
DRKAAKTLGIVMGVFLACWLPCFLAVLIDPYLDYSTPILDLLVWLRYF
NSTCNPLIHGFFNPWFQKAFKYIVSGKIFSSHSETANLFPEAH

Product Description

Expression Systems

Baculovirus/Insect expression system

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TAAR3P

Full Name

Trace amine associated receptor 3, pseudogene

Alternative Names

GPR57; GPR58; TAAR2; TAAR3; G protein-coupled receptor 57; TAAR3P; Trace amine associated receptor 3, pseudogene

Gene ID

[9288](#)

UniProt ID

[Q9P1P4](#)