

Product Information

MemDX™ Membrane Protein Human SERINC3 (Serine incorporator 3) Full Length

Cat. No.: **MPC0715K**

This product is for research use only and is not intended for diagnostic use.

This product is a 52.5 kDa Human SERINC3 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

SERINC3

Protein Length

Full length

Protein Class

Host-virus interaction; Immunity; Tumor target

Molecular Weight

52.5 kDa

TMD

8

Sequence

MGAVLGVFSLASWVPCLCSGASCLLCSCCPNSKNSTVTRLIYAFILLST
VVSyimQRKEMETYLKKIPGFCEGGFKIHEADINADKDCDVLVGYKAVYR
ISFAMAIFFVFVSLLMFKVKTSKDLRAAVHNGFWFFKIAALIGIMVGSFY
IPGGYFSSVWFVVGMI GAALFILIQLVLLVDFAHSWNESWVNRMEEGNPR
LWYAALLSFTSAFYILSIICVGLLYTYTKPDGCTENKFFISINLILCVV
ASIISHPKIQEHQPRSGLLQSSLITLYTMYLTWSAMSNEPDRSCNP NLM
SFITRITAPT LAPGNSTAVVPTPTPPSKSGSLLDSDNFIGLFVFLCLLY
SSIRTSTNSQVDKLTLSGSDSVILGDTTTS GASDEEDGQPRRAVDNEKEG
VQYSYSLFHLMLCLASLYIMMTLTSWYSPDAKFQSM TSKWPAVWVKISS
WVCLLLYVWTLVAPLVLT SRDFS

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

SERINC3

Full Name

Serine incorporator 3

Alternative Names

TDE; TDE1; AIGP1; TMS-1; DIFF33; SBBI99; serine incorporator 3; tumor differentially expressed protein 1; SERINC3; Serine incorporator 3

Gene ID

[10955](#)

UniProt ID

[Q13530](#)