

Product Information

MemDX™ Membrane Protein Human PLSCR3 (Phospholipid scramblase 3) Full Length

Cat. No.: **MPC4547K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human PLSCR3 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

PLSCR3

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MAGYLPPKGYAPSPPPYPVTPGYPEPALHPGPGQAPVPAQVPAPAPGFA
LFPSPGPVALGSAAPFLPLPGVPSGLEFLVQIDQILIHQKAERVETFLGW
ETCNRYELRSGAGQPLGQAAEESNCCARLCCGARRPLRVRLADPGDREVL
RLLRPLHCGCSCCPCGLQEMEVQAPP GTTIGHVLQ TWH PFLPKFSIQDAD
RQTVLRVVGPCWTCGCGTDTNFEVKTRDESRSVGRISKQWGGLVREALTD
ADDFGLQFPLDL DVRVKAVLLGATFLIDYMF FEKRGGAGPSAVTS

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

PLSCR3

Full Name

Phospholipid scramblase 3

Alternative Names

PLSCR3; PL scramblase 3; Ca^{2+} -dependent phospholipid scramblase 3; Phospholipid scramblase 3

Gene ID

[57048](#)

UniProt ID

[Q9NRY6](#)