

Product Information

MemDX™ Membrane Protein Human RHBDD3 (Rhomboid domain containing 3) Full Length

Cat. No.: **MPC4616K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human RHBDD3 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

RHBDD3

Protein Length

Full length

Protein Class

Receptor

TMD

5

Sequence

MHARGPHGQLSPALPLASSVLMMLMSTLWLVGAGPGLVLAPELLLDPWQV
HRLLTALGHTALPGLLLSLLLLPTVGWQQECHLGLRFLHASALLALAS
GLLAVLLAGLGLSSAAGSCGYMPVHLAMLAGEGHRPRRPRGALPPWLSPW
LLLALTPLLSSEPPFLQLLCGLLAGLAYAAGAFRWLEPSEERRLQVLQEGV
LCRTLAGCWPLRLLATPGSLAELPVTHPAGVRPPIPGPPYVASPDLSHW
EDSALPPPSLRPVQPTWEGSSEAGLDWAGASFSPGTPMWAALDEQMLQEG
IQASLLDGPAQEPQSAPWLSKSSVSSLRLQLERMGFTEQAVVALAATG
RVEGAVSLLVGGQVGTETLVTHGKGGAHSEGP

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

RHBDD3

Full Name

Rhomboid domain containing 3

Alternative Names

RHBDD3; PTAG; C22orf3; HS984G1A; rhomboid domain-containing protein 3; hypothetical protein HS984G1A; pituitary tumor apoptosis; Rhomboid domain containing 3

Gene ID

[25807](#)

UniProt ID

[Q9Y3P4](#)