

Product Information

TRITON X-100

α -[4-(1,1,3,3-Tetramethylbutyl)phenyl]- ω -hydroxy-poly(oxy-1,2-ethanediyl); TRITON X-100

Cat. No.: **MPD0170K**

This product is for research use only and is not intended for diagnostic use.

MPD0170K-1.png

This product is a detergent for membrane protein solubilization. It is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Type

Detergent

Classification

MB Reagents

Molecular Formula

t-Oct-C₆H₄-(OCH₂CH₂)XOH, x=9-10

CAS Number

9002-93-1

Formula Weight

647

CMC in H₂O

0.010 - 0.016%

Size

500 ml; 1 L

Storage

Aliquot and store at -20°C.