

# Product Information

## TWEEN 80

**Polyoxyethylene(80)sorbitan monolaurate; Poly(oxy-1,2-ethanediyl) derivs., (Z)-sorbitan mono-9-octadecenoate; TWEEN 80**

Cat. No.: **MPD0173K**

This product is for research use only and is not intended for diagnostic use.

MPD0173K-1.png

This product is a detergent for membrane protein solubilization. It is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Type

Detergent

### Classification

MB Reagents

### Molecular Formula

C<sub>64</sub>H<sub>124</sub>O<sub>26</sub>

### CAS Number

9005-65-6

### Formula Weight

1310

### CMC in H<sub>2</sub>O

1.6E-05

### Size

500 ml; 1 L

### Storage

Aliquot and store at -20°C.