

# Product Information

## TWEEN 40

**Polyoxyethylene sorbitan monolaurate; TWEEN 40**

Cat. No.: **MPD0174K**

This product is for research use only and is not intended for diagnostic use.

MPD0174K-1.png

This product is a detergent for membrane protein solubilization. It is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Type

Detergent

### Classification

MB Reagents

### Molecular Formula

C<sub>62</sub>H<sub>123</sub>O<sub>26</sub>

### CAS Number

9005-66-7

### Formula Weight

1284

### CMC in H<sub>2</sub>O

0.027 mM

### Storage

Aliquot and store at -20°C.