

Product Information

DLPE

1,2-dilauroyl -sn-glycero-3-phosphoethanolamine

Cat. No.: **MPL0014K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

Phosphadylethanolamine (PE)

Appearance

White powder

CAS Number

59752-57-7

Molecular Weight

579

Solubility

CHCl₃/MeOH/H₂O (65:35:8)

Purity

>99%

Size

500mg

Storage

-20°C

Storage Shelf Time

3 years