

Product Information

MemDX™ Membrane Protein Human TRIQK (Triple QxxK/R motif containing) Expressed *in vitro E.coli* expression system, Full Length

Cat. No.: **MPX1036K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human TRIQK membrane protein expressed *in vitro E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TRIQK

Protein Length

Full Length

Protein Class

Receptor

TMD

1

Sequence

MGRKDAATIKLPVDQYRKQIGKQDYKKTPILRATKLKAEAKKTAIGIKEVGLVLAAILALLLAFYAFFYLRLTTDVPDLDQDED

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

TRIQQ

Full Name

Triple QxxK/R motif containing

Alternative Names

TRIQQ; C8orf83; PRO0845; UPF0599; protein TRIQQ; triple repetitive-sequence of QXXK/R protein homolog; Triple QxxK/R motif containing

Gene ID

[286144](#)

UniProt ID

[Q629K1](#)