

# Product Information

## **MemDX™ Membrane Protein Mouse Triqk (Triple QxxK/R motif containing) Expressed *in vitro* E.coli expression system, Full Length**

Cat. No.: **MPX1038K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Triqk membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Mouse

#### Target Protein

Triqk

#### Protein Length

Full Length

#### Protein Class

Receptor

#### TMD

1

#### Sequence

MGRKDSSTKL PVDQYRKQIGKQDYKKT PILRATKLKAEAKKTAIGIKEVGLMLAAILALLLAFYAFFYLRLSTNIDSDLDLDED

### Product Description

#### Expression Systems

*in vitro* E.coli expression system

#### Tag

10xHis tag at the N-terminus

#### Protein Format

Soluble

#### Form

Liquid or Lyophilized powder

#### Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

### **Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

### **Target**

#### **Target Protein**

Triqk

#### **Full Name**

Triple QxxK/R motif containing

#### **Alternative Names**

Triqk; Tri; Gm11818; C130086A10; triple QxxK/R motif-containing protein; protein TRIQK; triple repetitive-sequence of QXXK/R protein homolog; Triple QxxK/R motif containing

#### **Gene ID**

[208820](#)

#### **UniProt ID**

[B2B9E1](#)