

Product Information

MemDX™ Membrane Protein Human OR4Q2 (Olfactory receptor family 4 subfamily Q member 2 (gene/pseudogene)) Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX3136K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human OR4Q2 membrane protein expressed *in vitro* *E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

OR4Q2

Protein Length

Full Length

Protein Class

GPCR

TMD

7

Sequence

MDKNQTEVMREFFLSGFSQTPSIEAGLFVLFLFFYMSIWVGNVLM/TVASDKYLNSSPMYFLLGNLSFLDLCYSTVTTPKLLADFFN

Product Description

Expression Systems

in vitro *E.coli* expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

OR4Q2

Full Name

Olfactory receptor family 4 subfamily Q member 2 (gene/pseudogene)

Introduction

Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. This olfactory receptor gene is a segregating pseudogene, where some individuals have an allele that encodes a functional olfactory receptor, while other individuals have an allele encoding a protein that is predicted to be non-functional.

Alternative Names

OR4Q2; OR4Q2P; olfactory receptor 4Q2; olfactory receptor OR14-21; Olfactory receptor family 4 subfamily Q member 2 (gene/pseudogene)

Gene ID

[390432](#)

UniProt ID

[P0C623](#)