

Product Information

MemDX™ Membrane Protein Human CMKLR2 (Chemerin chemokine-like receptor 2)

Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX3435K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human CMKLR2 membrane protein expressed *in vitro* *E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

CMKLR2

Protein Length

Full Length

Protein Class

GPCR

TMD

7

Sequence

MEDLEETLFEFENYSYDLDYYSLESDLEEKVQLGVVHWVSLVLYCLAFVLGIPGNAIVWFTGFKWKKTVTTLWFLNLAIDFIFLLF

Product Description

Expression Systems

in vitro *E.coli* expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

CMKLR2

Full Name

Chemerin chemokine-like receptor 2

Alternative Names

CMKLR2; GPR1; G-protein coupled receptor 1; probable G-protein coupled receptor 1; Chemerin chemokine-like receptor 2

Gene ID

[2825](#)

UniProt ID

[P46091](#)