

Product Information

MemDX™ Human MS4A1 A20 Cell Line

Cat. No.: **S01YF-0324-KX22**

This product is for research use only and is not intended for diagnostic use.

Product Information

Target Protein

MS4A1

Target Protein Species

Human

Host Cell Type

A20

Target Classification

Immune Checkpoint

Target Family

Immune Checkpoint

Target Research Area

Autoimmune Research; Immunology Research

Related Diseases

Immunodeficiency

Product Properties

Assay Types

Functional assay and biological assay

Assay Reporter

Luciferase

Mycoplasma Testing

Negative

Biosafety Level

Level 1

Activity

Expression verification; tumorigenicity verification

Form

Frozen cells

Freeze Medium

70% RPMI 1640 + 20% FBS + 10% DMSO

Culture Medium

RPMI 1640 + 10% FBS

Selective Antibiotic(s)

Regular antibiotics active against mycoplasmas, bacteria and fungi.

Handling Notes

Frozen cells should be thawed immediately upon receipt and grown according to handling procedure to ensure cell viability and proper assay performance.

Note: Do not freeze the cells upon receipt as it may result in irreversible damage to the cell line.

Disclaimer: We cannot guarantee cell viability if the cells are not thawed immediately upon receipt and grown according to handling procedure.

Incubation

37°C with 5% CO₂

Applications

Drug screening and biological assays

Application Notes

Cells were plated in a 384-well plate and incubated overnight at 37°C and 5% CO₂ to allow the cells to attach and grow. Cells were then stimulated with a control for high-throughput drugs screening and functional assays.

Use Restrictions

These cells are distributed for research use only.

Shipping

Dry ice

Storage

Liquid nitrogen

Target

Full Name

Membrane spanning 4-domains A1

Introduction

This gene encodes a member of the membrane-spanning 4A gene family. Members of this nascent protein family are characterized by common structural features and similar intron/exon splice boundaries and display unique expression patterns among hematopoietic cells and nonlymphoid tissues. This gene encodes a B-lymphocyte surface molecule which plays a role in the development and differentiation of B-cells into plasma cells. This family member is localized to 11q12, among a cluster of family members. Alternative splicing of this gene results in two transcript variants which encode the same protein.

Alternative Names

MS4A1; B1; S7; Bp35; CD20; FMC7; CVID5; MS4A2; LEU-16; B-lymphocyte antigen CD20; B-lymphocyte cell-surface antigen B1; CD20 antigen; CD20 receptor; leukocyte surface antigen Leu-16; membrane-spanning 4-domains, subfamily A, member 1; Membrane spanning 4-domains A1

Gene ID

931

UniProt ID

P11836