

Product Information

MemDX™ Membrane Protein Human CNR2 (1-33aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX134**

This product is for research use only and is not intended for diagnostic use.

This product is a Human CNR2 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

CNR2

Protein Length

Partial (1-33aa)

Molecular Weight

29.8 kDa

Sequence

1-33aa (Accession # P34972)

Product Description

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

CNR2

Full Name

Cannabinoid receptor 2

Alternative Names

CB2; CX5; CB-2; cannabinoid receptor 2 (macrophage); CNR2; Cannabinoid receptor 2

Gene ID

[1269](#)

UniProt ID

[P34972](#)