

Product Information

MemDX™ Membrane Protein Human IL13RA1 (23-343aa) Expressed in HEK293, C-6×His tag

Cat. No.: **S01YF-0423-KX149**

This product is for research use only and is not intended for diagnostic use.

This product is a Human IL13RA1 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

IL13RA1

Protein Length

Partial (23-343aa)

Molecular Weight

37.7 kDa

Sequence

23-343aa (Accession # P78552)

Product Description

Expression Systems

HEK293

Tag

C-6×His tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

IL13RA1

Full Name

Interleukin 13 receptor subunit alpha 1

Alternative Names

IL13RA1; NR4; CT19; CD213A1; IL-13Ra; interleukin-13 receptor subunit alpha-1; CD213a1 antigen; IL-13 receptor subunit alpha-1; IL-13R subunit alpha-1; IL13 receptor alpha-1 chain; cancer/testis antigen 19; interleukin 13 receptor, alpha 1; Interleukin 13 receptor subunit alpha 1

Gene ID

[3597](#)

UniProt ID

[P78552](#)