

Product Information

MemDX™ Membrane Protein Human IL5 (20-134aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX227**

This product is for research use only and is not intended for diagnostic use.

This product is a Human IL5 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

IL5

Protein Length

Partial (20-134aa)

Molecular Weight

39.3 kDa

Sequence

20-134aa (Accession # P05113)

Product Description

Activity

Yes

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

IL5

Full Name

Interleukin 5

Alternative Names

EDF; TRF; IL-5; interleukin-5; B-cell differentiation factor I; T-cell replacing factor; colony-stimulating factor, eosinophil; eosinophil differentiation factor

Gene ID

[3567](#)

UniProt ID

[P05113](#)