

Product Information

MemDX™ Membrane Protein Human PCSK9 (31-692aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX252**

This product is for research use only and is not intended for diagnostic use.

This product is a Human PCSK9 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

PCSK9

Protein Length

Partial (31-692aa)

Molecular Weight

97.2 kDa

Sequence

31-692aa (Accession # Q8NBP7)

Product Description

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

PCSK9

Full Name

Proprotein convertase subtilisin/kexin type 9

Alternative Names

FH3; PC9; FHCL3; NARC1; LDLQC1; NARC-1; HCHOLA3; proprotein convertase subtilisin/kexin type 9; convertase subtilisin/kexin type 9 preproprotein; neural apoptosis regulated convertase 1; subtilisin/kexin-like protease PC9

Gene ID

[255738](#)

UniProt ID

[Q8NBP7](#)