

Product Information

MemDX™ Membrane Protein Human GP6 (21-267aa) Expressed in HEK293, C-Fc tag

Cat. No.: **S01YF-0423-KX321**

This product is for research use only and is not intended for diagnostic use.

This product is a Human GP6 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

GP6

Protein Length

Partial (21-267aa)

Molecular Weight

53.1 kDa

Sequence

21-267aa (Accession # Q9HCN6)

Product Description

Expression Systems

HEK293

Tag

C-Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

GP6

Full Name

Glycoprotein VI platelet

Alternative Names

GP6; GPIV; GPVI; BDPLT11; platelet glycoprotein VI; glycoprotein 6; platelet collagen receptor; Glycoprotein VI platelet

Gene ID

[51206](#)

UniProt ID

[Q9HCN6](#)