

Product Information

MemDX™ Membrane Protein Human LGR4 (25-544aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX354**

This product is for research use only and is not intended for diagnostic use.

This product is a Human LGR4 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

LGR4

Protein Length

Partial (25-544aa)

Molecular Weight

83.3 kDa

Sequence

25-544aa (Accession # Q9BXB1)

Product Description

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

LGR4

Full Name

Leucine rich repeat containing G protein-coupled receptor 4

Alternative Names

GPR48; BNMD17; G protein-coupled receptor 48; LGR4; Leucine rich repeat containing G protein-coupled receptor 4

Gene ID

[55366](#)

UniProt ID

[Q9BXB1](#)