

Product Information

MemDX™ Membrane Protein Human GCGR (26-136aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX398**

This product is for research use only and is not intended for diagnostic use.

This product is a Human GCGR membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

GCGR

Protein Length

Partial (26-136aa)

Molecular Weight

39.2 kDa

Sequence

26-136aa (Accession # P47871)

Product Description

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

GCGR

Full Name

Glucagon receptor

Alternative Names

GGR; GL-R; glucagon receptor; GCGR; Glucagon receptor

Gene ID

[2642](#)

UniProt ID

[P47871](#)