

Product Information

MemDX™ Membrane Protein Human ENTPD1 (40-476aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX479**

This product is for research use only and is not intended for diagnostic use.

This product is a Human ENTPD1 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

ENTPD1

Protein Length

Partial (40-476aa)

Molecular Weight

75.9 kDa

Sequence

40-476aa (Accession # P49961)

Product Description

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

ENTPD1

Full Name

Ectonucleoside triphosphate diphosphohydrolase 1

Alternative Names

ENTPD1; CD39; SPG64; ATPDase; NTPDase-1; CD39 antigen; ecto-ATP diphosphohydrolase 1; ecto-ATPDase 1; ecto-ATPase 1; ecto-apyrase; lymphoid cell activation antigen; Ectonucleoside triphosphate diphosphohydrolase 1

Gene ID

[953](#)

UniProt ID

[P49961](#)