

Product Information

MemDX™ Membrane Protein Human IL17A (24-155aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX517**

This product is for research use only and is not intended for diagnostic use.

This product is a Human IL17A membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

IL17A

Protein Length

Partial (24-155aa)

Molecular Weight

41.3 kDa

Sequence

24-155aa (Accession # Q16552)

Product Description

Activity

Yes

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

IL17A

Full Name

Interleukin 17A

Alternative Names

IL17; CTLA8; IL-17; ILA17; CTLA-8; IL-17A; interleukin-17A; cytotoxic T-lymphocyte-associated antigen 8; cytotoxic T-lymphocyte-associated protein 8; interleukin 17 (cytotoxic T-lymphocyte-associated serine esterase 8); interleukin-17

Gene ID

[3605](#)

UniProt ID

[Q16552](#)