

Product Information

MemDX™ Membrane Protein Cynomolgus CXCL10 (22-98aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX805**

This product is for research use only and is not intended for diagnostic use.

This product is a Cynomolgus CXCL10 membrane protein expressed in HEK293.

Product Specifications

Host Species

Cynomolgus

Target Protein

CXCL10

Protein Length

Partial (22-98aa)

Molecular Weight

34.8 kDa

Sequence

22-98aa (Accession # F7H741)

Product Description

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

CXCL10

Full Name

C-X-C motif chemokine ligand 10

Alternative Names

C-X-C motif chemokine 10; CXCL10

Gene ID

[102134351](#)

UniProt ID

[A0A2K5WKN6](#)