

Product Information

MemDX™ Membrane Protein Mouse CD27 (24-182aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX826**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse CD27 membrane protein expressed in HEK293.

Product Specifications

Host Species

Mouse

Target Protein

CD27

Protein Length

Partial (24-182aa)

Molecular Weight

43.9 kDa

Sequence

24-182aa (Accession # P41272)

Product Description

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

CD27

Full Name

CD27 antigen

Alternative Names

S152; Tp55; Tnfrsf7; CD27 antigen; CD antigen 27; CD27L receptor; T-cell activation antigen CD27; tumor necrosis factor receptor superfamily member 7

Gene ID

[21940](#)

UniProt ID

[P41272](#)