

## Product Information

### **MemDX™ Membrane Protein Mouse TNFRSF8 (19-287aa) Expressed in HEK293, C-Human Fc tag**

Cat. No.: **S01YF-0423-KX834**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse TNFRSF8 membrane protein expressed in HEK293.

#### Product Specifications

##### Host Species

Mouse

##### Target Protein

TNFRSF8

##### Protein Length

Partial (19-287aa)

##### Molecular Weight

54.4 kDa

##### Sequence

19-287aa (Accession # Q60846)

#### Product Description

##### Expression Systems

HEK293

##### Tag

C-Human Fc tag

##### Protein Format

Soluble

##### Form

LYOPH

##### Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

##### Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

**Storage**

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

**Target****Target Protein**

TNFRSF8

**Full Name**

Tumor necrosis factor receptor superfamily, member 8

**Alternative Names**

Ki; Cd30; Ki-1; D1S166E; tumor necrosis factor receptor superfamily member 8; CD30L receptor; lymphocyte activation antigen CD30

**Gene ID**

[21941](#)

**UniProt ID**

[Q60846](#)