

Product Information

MemDX™ Membrane Protein Mouse HBEGF (24-160aa) Expressed in HEK293, C-6×His tag

Cat. No.: **S01YF-0423-KX835**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse HBEGF membrane protein expressed in HEK293.

Product Specifications

Host Species

Mouse

Target Protein

HBEGF

Protein Length

Partial (24-160aa)

Molecular Weight

16.0 kDa

Sequence

24-160aa (Accession # Q06186)

Product Description

Expression Systems

HEK293

Tag

C-6×His tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

HBEGF

Full Name

Heparin-binding EGF-like growth factor

Alternative Names

Dtr; Dts; Hegfl; proheparin-binding EGF-like growth factor; diphtheria toxin receptor; heparin binding epidermal growth factor-like growth factor

Gene ID

[15200](#)

UniProt ID

[Q06186](#)