

Product Information

MemDX™ Membrane Protein Mouse MSLN (298-600aa) Expressed in HEK293, C-6×His tag

Cat. No.: **S01YF-0423-KX855**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse MSLN membrane protein expressed in HEK293.

Product Specifications

Host Species

Mouse

Target Protein

MSLN

Protein Length

Partial (298-600aa)

Molecular Weight

35.0 kDa

Sequence

298-600aa (Accession # Q61468)

Product Description

Expression Systems

HEK293

Tag

C-6×His tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

MSLN

Full Name

Mesothelin

Alternative Names

MPF; mesothelin; megakaryocyte potentiating factor; pre-pro-megakaryocyte-potentiating factor

Gene ID

[56047](#)

UniProt ID

[Q61468](#)